

FEATURES

South-Tek's Ultra Zero Air Generator Series produce dry and hydrocarbon / NOx-free air, using air from oil free compressors, thus avoiding the need to use traditional bottles that are often complex to change.

- | | |
|--|---|
| <ul style="list-style-type: none"> Increased detector sensitivity ensures precise, analytical results Easy Installation; connect a suitable oil free compressor (optional) to almost immediately have air at the right grade Continuous operation for improved lab efficiency Can be installed on lab wall | <p>Ideal for:</p> <ul style="list-style-type: none"> GC-FID GC-NPD GC-FPD GC-PFPD THA TOC |
|--|---|



UNIT DETAILS

	GT 1500 plus	GT 3000 plus	GT 6000 plus	GT 15000 plus	GT 30000 plus
Dimensions					
A (Width)	21.6"			24.8"	
B (Depth)	9"			12.2"	
C (Height)	16"			18.5"	
Weight	22 lbs.	33 lbs.		56 lbs.	
Unit Details					
Flow l / min.	1.5	3	6	15	30
NOx, SOx, O ₃ , HC & CO out	< 0.1				
CO ₂ Out	< 5 ppm				
Max CO in	50 ppm				
Max HC in	100 ppm				
Dew point	50°C				
Pressure in (*)	4.5-10 bar				
Pressure drop	< 1 bar				
Stability (min)	45				
In / out connections	1 / 4-1 / 8 (NPT)				
Power	110-120V 60Hz / 220-240V 50Hz				
Working temp	Amb.+15°C - 35°C				
Warranty	2 Year Warranty ¹				

¹. Per Manufacturer's Terms & Conditions.