

FEATURES

South-Tek's Zero Air Generator Series designed with safety and efficiency in mind, produces dry and hydrocarbon-free air, using air from oil-free compressors.

- | Increased detector sensitivity ensures precise, analytical results
- | Ergonomic, compact and silent design
- | Easy Installation
- | Continuous operation for improved lab efficiency
- | European plugs available upon request

- Ideal for:
- | GC-FID
 - | GC-NPD
 - | GC-FPD
 - | GC-PFPD
 - | THA



UNIT DETAILS

	GC 1500 plus	GC 3000 plus	GC 6000 plus	GC 15000 plus	GC 30000 plus
Dimensions					
A (Width)	21.6"				
B (Depth)	9"				
C (Height)	16"				
Weight	21 lbs.	24 lbs.		26.5 lbs.	28.5 lbs.
Unit Details					
Flow l / min.	1.5	3	6	15	30
HC & CO out	< 0.1				
Max CO in	50 ppm				
Max HC in	100 ppm				
Pressure in (*)	4.5-10 bar				
Pressure drop	< 1 bar				
Stability (min)	45				
In / out connections	1 / 4-1 / 8 (NPT)				
Power	110-120V 60Hz / 220-240V 50Hz				
Working temp	Amb.+15°C - 35°C				
Warranty	2 Year Warranty ¹				

1. Per Manufacturer's Terms & Conditions.